

HPP Observatory Quick User Guide

August 2015

Abstract

This step-by-step user guide describes the standard procedure for start-up, operation, and shutdown of the HPP observatory as it is adequate for most student projects (e.g. VP-ASL experiments) or for outreach activities. For more functionalities of the telescope, the mount, or the camera see the corresponding manuals.

Contents

1	Instrument specifications	3
2	Start-up	4
2.1	Open the telescope box (“the dome”)	4
2.2	Remove mirror cover	4
2.3	Mounting the SBIG camera	4
2.4	Power-on the instruments	4
2.5	Turn-on GM4000 mount	5
2.6	Start-up the instrument software	5
2.7	Figures	6
3	Shutdown	17
4	Instrument operation using individual software packages	18
4.1	GM4000 mount operation	18
4.1.1	Alignment of the telescope	19
4.2	CDK operation	19
4.3	Camera operation (SBIG STX-16803)	19
4.4	Figures	21
5	Instrument operation using MaxImDL	25
5.1	Observatory control	25
5.1.1	Setup	25
5.1.2	Telescope control	25
5.1.3	Focusing options	26
5.2	Camera control	27
5.3	Image inspection	28
5.4	Figures	28
6	Spectroscopy with the DADOS spectrograph	36
6.1	General comments	36
6.2	Spectrograph assembly and setup	36
6.3	Alignment	37
6.3.1	Spatial orientation	37

6.3.2	Focal alignment	37
6.4	Operation	38
6.5	Calibration (TBD)	38
6.6	Figures	38
7	Data reduction using MaxImDL (TBD)	47
A	Mounting the SBIG camera	48
B	First light of DADOS spectrograph	53

1 Instrument specifications

telescope <i>Planewave CDK20</i>	
optical design	Corrected Dall-Kirkham
aperture	20 inch = 508 mm
focal length	3454 mm
diffraction limit at $\lambda = 500$ nm	0.2 arcsec
focal ratio	F/6.8
central obstruction	39 % of the primary mirror diameter
back focus from racked in focuser	147 mm
camera <i>SBIG STX-16803</i>	
CCD	KAF-16803
array	4096×4096 pixels ²
pixel size	$9 \times 9 \mu\text{m}^2$
CCD dimensions	36.9×36.9 mm ²
CCD type	full frame, antiblooming
shutter	even-illumination mechanical shutter
external ports	USB20/1.x and 10/100 Mbs Ethernet
telescope + camera	
field of view	36.7×36.7 arcmin ²
diagonal field of view	51.9 arcmin
plate scale	0.54 arcsec/pixel
magnitude limit	TBC

2 Start-up

Make sure that you go through all of the following operations, and that you stick to the order.

2.1 Open the telescope box (“the dome”)

In general, observations are only allowed if the dome can be completely opened and locked in its open position. In the winter make sure that no snow or dirt can fall down onto the telescope when opening/sliding the dome! Only push/pull the dome from the inside so that you can check that the telescope fits through the door opening!

- fully open both doors
- make sure the open doors are snapped into their secure locks (they cannot swing back, see Fig. 1)
- release the dome locks (located inside the dome, see Fig. 1)
- check that the telescope is correctly parked and that it will fit through the doors. (if this is not the case, stop all activities and call the HPP responsible or the VP-ASL supervisor at IfA)
- push the dome away into its open position
- engage the dome locks
- close the doors

2.2 Remove mirror cover

When the telescope is parked, in particular when it is not used for a long time period, the primary mirror should be protected by the corresponding plexiglas cover (see Fig. 2). For removing/installing the cover pull back the spandex stray light cover and pass the mirror cover through the bars as shown in Fig. 2.

2.3 Mounting the SBIG camera

The SBIG camera should always be attached to the telescope. If this is not the case or if you intent to use the telescope without camera, see Appendix A for instructions about how to mount or dismount the camera.

Important: if the camera is detached use the strap at the telescope handle to fix the cables! (see Fig. 3)

2.4 Power-on the instruments

In general, all required instruments except the remote computer will be powered-on by activating the main power switch (see Figs. 4, 5).

Default setting:

- remote computer - OFF
- power adapter GM4000 mount - ON
- Delta-T control box - ON
- Delta-T heater switches (Fig. 8, primary, secondary) - OFF (LEDs off)
- CDK Power (Fig. 9, 10) - ON
- CDK fan switch (Fig. 10)- OFF
- SBIG - ON

The power switch of the remote computer is on the backside (see Fig. 6).

Note that the SBIG does not have an on/off switch. It will start-up and initialize as soon as it is connected to the power supply. Immediately after activating the main power switch you should hear the fan of the camera.

2.5 Turn-on GM4000 mount

The GM4000 is a computerized mount and needs to be booted (and shut down after use! see Sect. 3). To start the mount press the on/off button at the GM4000 control box (Figs. 4, 7). Note that although the button looks like a switch it is actually only a button and it will resume its initial position after pressing it (so press the button only once!). After a few seconds the display of the GM4000 hand control will consecutively read “booting”- “starting-up”- “version”- “initializing”. Initialization takes about 30 seconds and the mount is ready to use as soon as the display reads the time and date.

Note that for visual observations the telescope is now ready and can be operated by the hand control. For an overview of the basic functionalities of the hand control see Sect. 4.1.

2.6 Start-up the instrument software

All instruments can be fully operated by software installed on the remote computer outside at the telescope. The computer can be remotely accessed by the control room computer. Turn on both computers. To connect to the remote computer click on the *Remote Desktop* shortcut (see Fig. 11) and use the login provided by the ASL supervisor. Ignore the warning after clicking on *Connect* and click on *Yes*.

- control room:
username: XXXXXXXX
password: XXXXXXXX
- telescope computer:
username: XXXXXXXX
password: XXXXXXXX

Important: make sure to shutdown the remote computer before turning off the main switch (see Sect. 3 for the correct shutdown procedure!)

You find the relevant programs grouped together on the desktop of the remote computer (see Fig. 12):

- main program: **MaxImDL**
this software is used to control all individual components
- auxiliary software:
these are the individual software distributions from the manufacturers but they are obsolete when using MaxImDL.
mount control: **Virtual Keypad**
telescope control / focus control: **PWI3**
camera control: **CCDOps**

See the next sections 4 and 5 for operational instructions of the software.

2.7 Figures




Figure 1: door locks (*top*) and dome locks (*bottom*).




Figure 2: Mirror cover.



Figure 3: Strap to fix the loose cable when camera is detached.



Figure 4: Overview of electronics: main power switch (circle), remote computer (A), power adapter GM4000 mount (B), router (B), GM4000 control box (C), Delta-T control box (D), (handles to push/pull dome (E)).




Figure 5: Main power switch (circle). Also indicated are the remote computer (A) and the power adapter for the GM4000 mount (B).



Figure 6: Power switch of the remote computer (circle).




Figure 7: On/Off button at the GM4000 control box.




Figure 8: Delta-T control box. Frontside (*top*) and backside (*bottom*). The circles indicate the heater switches in their default positions (OFF, LEDs off).



Figure 9: CDK electronics overview with the CDK control panel on top and the Delta-T control box at the bottom (see also Fig. 8).




Figure 10: CDK control panel with default setting: Power ON (A), Fan OFF (B).




Figure 11: Shortcut and login for remote access to telescope computer. Note: ignore the warning after clicking *Connect* and click on *Yes*.



Figure 12: Software shortcuts on remote computer.




Figure 13: Finder scope. Off position (circle) for the light modulator (A) from the front (top) and the back (bottom). Also indicated are the caps (B) which protect the alignment knobs for x and y direction.

3 Shutdown

!!! To prevent any damage to the instrument and the software it is very important to follow a correct shutdown procedure !!!

For shutting down reverse the start-up procedure:

- **Shutdown Camera** Disconnect the camera. In **MaxImDL** go to the camera control window and click on the *Warm up* button. Wait until duty cycle is 0 % then click on the *Disconnect* button. In **CCDOps** use the *ShtDn* quick launch button. In the pop-up panel then select *Warm Up & Wait*. After a few seconds the software will disconnect the link with the camera and you can close CCDOps.
- **Disconnect observatory utilities** In **MaxImDL** go to the observatory control window and individually disconnect the **webcam**, the **focuser**, and the **telescope**. Then close the **PWI3** then close the programs. If the **Virtual keypad** was used disconnect it and close the program.
- Turn off the light of the **Finder scope** (see Fig. 13)
- **Park the GM4000 mount** (see also comments in Sect. 4.1)
- **Shutdown the GM4000 mount** press the on/off button at the GM4000 control box (Figs. 4, 7). Note that although the button looks like a switch it is actually only a button and it will resume its initial position after pressing it (**so press the button only once!**). **After up to 5 seconds** the mount will start the shutdown procedure and the display will read *shutting down*. Wait until the display gets dark.
- **Shutdown the remote computer** by **alt+F4** and select *Herunterfahren*
- **turn off main power switch** (Figs. 4, 5)
- install mirror cover (see Sect. 2.2, Fig. 2)
- open the dome doors and **lock them in the secure locks**
- release the dome locks
- close the dome **from the inside** so that you can check that the telescope fits through the door opening (i.e. that the correct parking position was set).
- **engage the dome locks**
- close doors and lock the dome

4 Instrument operation using individual software packages

This section gives a short overview of basic functionalities of the instrument software. These should be sufficient for most applications. For more specialized instrument operations see the original user manuals of the instruments.

4.1 GM4000 mount operation

The GM4000 mount is operated by the **Virtual Keypad** software which is effectively only a virtualization of the physical GM4000 hand control. Both the software and the hand control have exactly the same functionalities and menu navigations. After connection of the software (menu *Connection-Connect*) the physical hand control and the software are synchronized and therefore one can still use both (but not at the same time). For a complete overview of the functionalities of the GM4000 mount see the manual, in particular Chapt. 13 (*Menüstruktur*) which is also given in the document *GM4000_menu-structure*.

general

In general navigation through the hand control menus is very intuitive. Start by hitting *MENU* and then use the up/down keys (not the *E-W/N-S* keys!) to navigate through the menu, use *ENTER* for selecting submenus and for initializing actions, and use *ESC* for returning to the superior (sub-)menu. Starting from the main screen (outside the menu, the display reads the RA and DEC coordinates) you can use the quick access/action buttons as indicated on the keys to directly enter various object catalogs for selecting new targets and for using the automatic goto options. The *E-W/N-S* keys can be used for manually slewing the telescope to a given position (after unparking the mount, see below). To change the slew speed go to the main screen and use the up/down keys.

Note when using the automatic goto function the mount will indicate by a beep signal when it has arrived the new position (turn on the PC-speakers when you are using the virtual keypad). The fine positioning can sometimes take a few seconds therefore wait with new operations until you hear the signal.

Note that when the telescope is moving automatically (i.e. by using a goto function) it will immediately stop as soon as you hit one of the *E-W/N-S* keys or the *STOP* key.

unparking / parking of the telescope

For storage inside the dome the telescope needs to be in a fixed parking position predefined in the hand control. As long as the mount is in the parked state it is not possible to move the telescope.

- press *MENU*
- select *ALIGNMENT* (by using the up/down keys, not the *N-S/E-W* keys!) and hit *ENTER*
- select *UNPARK/PARK* and hit *ENTER* twice

Note that in case of unparking the mount will simply change its state but the telescope will stay in its parking position. In case of parking the telescope will automatically move to its parking position (the parking procedure can be started from any telescope orientation).

tracking speed

Depending on the object you need to change the tracking speed:

- *MENU-DRIVE-TRACKING SPEED*

One can choose between *sidereal*, *solar* (for solar system, e.g. planets), *lunar*, *custom* (only for advanced users), and *stop* (no tracking).

4.1.1 Alignment of the telescope

If the pointing of the telescope is bad you can either try to load a different sky model

- *MENU-ALIGNMENT-ALIGN DATABASE-LOAD MODEL*

or you have to redo the alignment procedure:

- *MENU-ALIGNMENT-3-STARS* (follow the instructions on the display)

For the alignment procedure the GM4000 mount requires calibration/centering of 3 reference stars plus a few additional stars scattered on the sky for refining the sky model. From the list of available stars (suggested by the software) choose bright stars that you know and that are preferably widely scattered on the sky. The telescope will then slew to the (offset) coordinates of the selected star and you will need to use the *N-S/E-W* keys to center the star on the CCD. After centering the star confirm it and continue with the next star.

Note that for the first 3 reference stars it is not possible to cancel the calibration and to choose a different star if the initially chosen star is not visible (e.g. because of clouds) or not found. After calibration of the 3 reference stars add additional stars for refining the correction model. Here you can abort if the star is not visible. Usually about 6 stars (3 reference stars + 3 additional stars) are sufficient enough for accurate pointing. For later use do not forget to save the new model by

- *MENU-ALIGNMENT-ALIGN DATABASE-SAVE MODEL.*

4.2 CDK operation

The software **PWI3** is used for setting/monitoring the mirror temperature and the mirror fan, as well as for operating the automated focusing stage. For the moment the software should only be used for focusing. After launching the software click on the button *CONNECT* as indicated in Fig. 14. If the connection was successful the software will immediately start monitoring and displaying the ambient and mirror temperatures.

Mirror temperature control

Not yet tested, do not use it for the moment.

Focusing

In the right side panel select the tab *Focuser* (Fig. 14) to operate the focusing stage. The actual position is displayed in μm . You can change the position by using the *IN/OUT* buttons or you can manually set another absolute position and hit *GOTO*. Note when setting a new position the motor will slightly overshoot and then slowly settle at the given position.

SBIG Filter	focus position at T=?TBD? [μm]
<i>R</i>	TBD
<i>G</i>	TBD
<i>B</i>	TBD
<i>L</i>	TBD
<i>C</i>	TBD
<i>Hα</i>	TBD

4.3 Camera operation (SBIG STX-16803)

The camera is operated by the software **CCDOps**. After launching the software the camera first needs to be linked to the software. Either use *Camera-Establish COM Link* or go directly to the menu *Camera-Setup* which will both link the camera and open a panel for setting the camera configuration (e.g. CCD temperature, binning mode, readout mode). Both options are also available from the quick launch tab shown in Fig. 15. When linking of the camera was successful this will be indicated in the monitoring panel at the bottom right together with the actual CCD temperature, the duty cycle to stabilize the temperature, and the selected filter (see Fig. 15).

standard operation settings (*Camera-Setup*, see *Fig. 16*)

- **CCD cooling** In order to reduce dark current the CCD can be cooled by a Peltier-element. Depending on ambient temperature use a setpoint between -20°C and 0°C . As soon as the temperature is reached the duty cycle should drop from 100 % to a value < 75 %. If this is not the case then the camera cannot stabilize the temperature and one should choose a higher setpoint.
- **resolution mode** This sets the on-chip binning value: *high* for no binning, *medium* for 2×2 binning, and *low* for 3×3 binning. Besides for increasing the flux (but lowering resolution) and therefore minimizing the integration time, binning is also used for minimizing the overhead time for downloading the image. In *high* resolution mode the image size is 4096×4096 pixels with a download time of TBC whereas for *low* resolution the image size is 1365×1365 pixels with a download time of TBC.
- **other settings** do not change other settings but use default values given in *Fig. 16*.

selection of filters

For selecting different filters go to the menu *Filter* and click on the new filter. The monitoring panel at the bottom right will display *Filter:moving* until the new filter is set.

acquiring images

For acquiring images use the quick launch button **Grab** (see *Fig. 15*). In the pop-up panel one can then choose the following options:

- **Exposure Time** Sets the integration time. Note that the shortest possible integration time is 0.12 s. When choosing an integration time < 0.12 s the software will automatically reset it to 0.12 s. Also note that for high precision photometry one should use integration times > 0.5 s or even > 1 s because the shutter might not be reliable for shorter integration times.
- **Dark Frame** Choose *Only* for recording of dark images only, *None* for recording of light images only, and *Also* for recording both light and dark images.
- **Image Size** Sets the (sub)image to be read out. Default is *Full* for reading out the complete CCD-array.
- **Exposure delay** not used, choose 0.
- **Special Processing** Choose *None* for single image recording and use *Auto Grab* for recording sequences of images. For single recording one can save the image after acquisition by using the quick launch button *Save* (see *Fig. 15*). If *Auto Grab* is selected a new pop-up window will appear where one can select the root filename and directory using the button *Set Name/Dir*, and the *Number of exposures* to be consecutively acquired. Do not change the other settings (use default values given in *Fig. 17*).

focusing mode

For focusing one can use the **Focus** button (see *Fig. 15*) for consecutively acquiring images. In the pop-up panel one can then choose the following options:

- **Exposure Time** Sets the integration time
- **Frame size** (Sub)array to be readout (e.g. to minimize overhead/download time)
- **Turbo Mode** Fast readout of CCD without proper resetting of the charges after readouts (e.g. to minimize overhead time, default is *Off*).
- **other settings** do not change other settings but use default values given in *Fig. 18*.

image inspection

For quick image inspection use the quick launch buttons *XHair* (cross-hairs), *Hist* (histogram), and/or use the *Contrast* pop-up panel which appears after each readout (see Fig. 15). In the *Contrast* panel one can manually set the greyscale range and select different zoom options (*Mag*).

4.4 Figures




Figure 14: PWI3 panel. Only use functions indicated by circles.




Figure 15: CCDOps main panel. Quick launch options and monitoring panel at the bottom right are indicated by circles. Open sub-panels: *cross-hairs*, *histogram*, and *contrast*.




Figure 16: Camera setup panel with default values and user options indicated by circles.




Figure 17: Autograb panel with default values and user options indicated by circles.




Figure 18: Focus panel with default values and user options indicated by circles.

5 Instrument operation using MaxImDL

The observatory software MaxImDL is based on the widely used ASCOM driver standards for telescopes, cameras, webcams, focuser assemblies, and many other instruments used in Astronomy. Therefore, all the individual components can be controlled by MaxImDL and advanced options, e.g. automatic focusing or planetarium utilities, can be used. In addition the software has a wide range of built-in data reduction utilities and can manage calibration data flows based on fits-header informations.

This section gives a short overview of basic functionalities of MaxImDL. These should be sufficient for most applications. For a first glance into the functionalities of MaxImDL see also the video tutorials (http://www.cyanogen.com/maxim_tut.php) provided by the manufacturer and/or the MaxImDL user manual.

Note: operation using MaxImDL assumes a well aligned telescope. If the pointing is bad or insufficient repeat the alignment procedure using the hand control as described in Section. 4.1.1.

5.1 Observatory control

All observatory utilities except the operation of the CCD camera (e.g. telescope pointing, focusing, planetarium utilities; for the camera operation see Sect. 5.2) are controlled by the **Observatory** window (e.g. see Figs. 20-27). The observatory window can be launched by the quick access button in the main panel of MaxImDL (see Fig. 19).

5.1.1 Setup

Go to the **Setup** tab (see Fig. 20). Connect each component separately starting from the top using the corresponding **Connect** buttons. For disconnecting start from the bottom.

Note:

When connecting the *Focuser* the **PWI3** software is launched in the background. Leave it open and only close it after closing MaxImDL.

Note:

The webcam is only used when using the Spectrograph.

5.1.2 Telescope control

The telescope can be controlled manually by the **Telescope** tab (see Fig. 21) or one can use the **Go to** button for known objects given by the **Catalog** tab (see Fig. 23) or by using the right-click option of the mouse by the **All Sky** tab (see Fig. 24) or the **Zoom** tab (see Fig. 25).

Telescope

The information window summarizes the current position, the tracking status and if the telescope is moving or Idle.

Use the *Mount* panel to **Park/Unpark** the telescope and to set **sidereal tracking** on/off. It is only possible to enable/disable *Sidereal Tracking / No Tracking* but setting the *Solar Tracking / Lunar Tracking* is only possible using the hand control (see Section 4.1).

The telescope can be moved manually in a similar way as with the keypad on the hand control by using the keys in the **Nudge** panel. The speed (the movement corresponding to one mouse click) can also be set there.

To move the telescope to specific coordinates use the **Target Coordinates** panel and click on the **Go To** button to move the telescope.

All the other options in this panel should already be correctly set. For example the **Site and Optics** settings are directly loaded from the GM4000 mount (see Fig. 22). Note that *Alt/Az* can only be displayed by the different panels if the *Site and Optics* settings are correctly set.

Catalog

To find a specific object one can either display objects using **Search ID** for searching for the object's name or one can restrict the list of all available objects by defining different Criteria using the **Search Region** option. For example by specifying a minimum altitude one can only display all currently visible objects.

To move the telescope to an object highlight the object in the list and click on the **Go to** button.

All Sky and Zoom

This tab shows a map of the sky. Use the **Options** button to change the *density* or *type* of the displayed objects and turn on/off the *object labels*. See also the mouse right-click and scrolling options, e.g. to zoom in/out or to move the telescope to the current mouse position.

5.1.3 Focusing options

The automated focuser is operated by the **Focus** tab in the observatory window (see Fig. 26). The **Focuser Status** panel displays the current focuser position and the results of the last focus measurement.

There are two approaches for focusing the image onto the camera:

Manual focusing

Use both the **Camera Control** window (see Sect. 5.2) and the **Observatory/Focuser** window:

- Take multiple exposures (see Sect. 5.2 and Fig. 31) while moving the focuser in/out by using either the **Incremental** or the **Absolute** panel (see Fig. 26).
To speed up the image acquisition choose a subframe around the object to be focussed (see Sect. 5.2).
- Optimize the *FWHM* or the *Half Flux Dia.* values given by the information panel in the Camera Control window.

Automatic focusing

Use the **Autofocus** panel in the **Observatory** window (see Figs. 26, 27):

- click on the **Options** button and set the options as given by Fig. 26.
- click on the **Snapshot** button to acquire a full-field image (using the current settings given in the **Camera Control** window).
- in the full-field snapshot image MaxImDL automatically selects the brightest star as focus reference. If another star should be used for focusing click on this star. The Information panel in the **Observatory/Focus** window will then display the coordinates of the chosen reference star.
- click on the **Exposure** button to define the exposure settings for automatic focusing (see Fig. 27), i.e. define the *Exposure Time*, the *Binning* mode, the size of the *Subframe*, and the *Filter*.
- Finally, click on the **Start** button in the **Autofocus** panel to start the automatic focusing procedure.
- MaxImDL then automatically acquires images of the chosen reference star for different focus positions. The results are continuously displayed in the Information panel and a *HFD / Position* plot is created (see Fig. 26).
- As soon as the minimum of the *HFD / Position* plot can be determined click on the **Abort** button to stop the focusing.
- Use the **Absolute** panel to move the focuser to the Minimum position.

5.2 Camera control

All camera operations (e.g. *setup* and *image acquisition*) are controlled by the **Camera Control** window (e.g. see Figs. 28, 31). The camera control window can be launched by the quick access button in the main panel of MaxImDL (see Fig. 19).

Setup

Go to the **Setup** tab (see Fig. 28). Click on **Setup Camera** to set the camera model. Use the settings according to Fig. 28 and only change the model (i.e. STX-16803 for standard imaging or ST-1603 for spectroscopy (Sect. 6)). Click on **Setup Filter** to set the filter settings (e.g. filter names, see Fig. 30). When using the ST-1603 camera set *Filter or Controlling Camera Model* to *No Filter* else use *SBIG Universal*. Note that the filter settings can only be accessed and changed when the camera is disconnected. If the filter settings are still ok (default situation) click on the **Connect** button. The camera should now be connected and some basic information (e.g. cooling status) is displayed in the information area (see Fig. 28). Activate the cooling by clicking on **Coolers on** and set a setpoint value by clicking on **Cooler** for Camera 1.

- **CCD cooling** In order to reduce dark current the CCD can be cooled by a Peltier-element. Depending on ambient temperature use a setpoint between -20°C and 0°C . As soon as the temperature is reached the duty cycle should drop from 100 % to a value $< 75\%$. If this is not the case then the camera cannot stabilize the temperature and one should choose a higher setpoint.

Image acquisition

For acquiring images go to the **Expose** tab (see Fig. 31). Here one can define all exposure settings and also define named presets for often used settings.

In particular use the following options:

- **Exposure Preset** Sets basic settings for commonly used camera presets (e.g. *Single Imaging*, *Continuous Imaging*, *Dark*). Choose the preset that matches your intention and adapt some of the settings below to optimize the exposure.
Note: whenever one of the pre-defined settings of the chosen preset is changed this is indicated by a star (e.g. **Single Imaging*).
- **Seconds** Sets the integration time. Note that the shortest possible integration time is 0.12 s. When choosing an integration time < 0.12 s the software will automatically reset it to 0.12 s. Also note that for high precision photometry one should use integration times > 0.5 s or even > 1 s because the shutter might not be reliable for shorter integration times.
- **Frame Type** Choose *Light* for recording normal images with shutter open and *Dark* for dark images with the shutter closed.
- **Filter Wheel** Select different filters. The monitoring panel at the bottom right will display *Waiting for Filter wheel* until the new filter is set.
- **Binning** Set the on-chip binning value (1 for no binning). Besides for increasing the flux (but lowering resolution) and therefore minimizing the integration time, binning is also used for minimizing the overhead time for downloading the image.
- **Subframe** Sets the (sub)image to be read out. Click *On* and draw the subframe using the mouse.
Note: for the ST-1603 Subframe does not work and needs to be de-activated!
- **Acquisition settings** Choose *Single* for single exposures, *Continuous* for multi-exposures, or *Autosave* for multiple exposure settings (see also Fig. 32).
- more **Options** (e.g. automatic dark exposure/subtraction) are available using the *Options* button.
- **Start / Stop** Start or stop the exposure with the chosen settings.

5.3 Image inspection

For quick image inspection use the stretch button (see Fig. 19) to manually set the greyscale range and use the **mouse-click options**:

- **scrolling** for zooming in/out
- **dragging + Shift + left-click** for adjusting brightness and contrast (similar to DS9 functionality)
- **dragging + Ctrl + left-click** for shifting the image

For more options see Sect. 7 (TBD).

5.4 Figures




Figure 19: Quick launch buttons for camera control and observatory control (two buttons on the right) and for setting the image cuts (i.e. *stretch*, left button.)




Figure 20: **Observatory:** Setup tab. Connect/Disconnect each component separately. The Webcam is only used for Spectrometry.




Figure 21: **Observatory:** Telescope tab.






Figure 22: **Observatory:** Telescope tab with *Site and Optics* sub-panel.

Observatory

All Sky | **Catalog** | Telescope | Dome | Focus | Status | Weather | Webcam | Setup

☐ Search ID ☒ Search Region Keypad  Options 

Object ID Search

Category:

Object ID:

Region Search

Criterion	Min	Max
Altitude	30.0	90.0

Search Results: found 8114 objects

ID1	ID2	Type	RA	Dec	Mag	Size	HA	Alt	Az	SlewDist
Albireo	Bet1 Cyg	Star	19h 30m 43s	27° 57' 34"	3.2	20h ...	45.4°	96.6°	130.8°	
Alderamin	Alp Cep	Star	21h 18m 35s	62° 35' 08"	2.4	18h ...	45.4°	40.0°	125.2°	
Alkaid	Eta UMa	Star	13h 47m 32s	49° 18' 48"	1.9	02h ...	67.8°	287.2°	71.4°	
Alphecca	Alp CrB	Star	15h 34m 41s	26° 42' 53"	2.2	00h ...	68.6°	196.8°	78.9°	
Arcturus	Alp Boo	Star	14h 15m 40s	19° 10' 57"	0.0	01h ...	54.4°	226.5°	59.2°	
Cor Caroli	Alp2 CVn	Star	12h 56m 02s	38° 19' 05"	2.9	03h ...	55.3°	272.0°	57.1°	
Deneb	Alp Cyg	Star	20h 41m 26s	45° 16' 49"	1.3	19h ...	43.4°	65.3°	135.8°	
Dubhe	Alp UMa	Star	11h 03m 44s	61° 45' 05"	1.9	04h ...	47.2°	317.7°	71.9°	
Etamin	Gam Dra	Star	17h 56m 36s	51° 29' 20"	2.2	22h ...	71.1°	66.9°	108.7°	
Kochab	Bet Umi	Star	14h 50m 42s	74° 09' 20"	2.1	01h ...	62.2°	349.6°	92.2°	
Kornephoros	Bet Her	Star	16h 30m 13s	21° 29' 23"	2.8	23h ...	63.4°	165.1°	89.9°	
Merak	Bet UMa	Star	11h 01m 50s	56° 22' 55"	2.4	04h ...	45.2°	310.4°	66.5°	
Mizar	Zet UMa	Star	13h 23m 56s	54° 55' 31"	2.2	02h ...	64.5°	301.7°	72.9°	
Polaris	Alp UMi	Star	02h 31m 46s	89° 15' 52"	2.1	13h ...	46.8°	0.3°	100.5°	
Ras Algethi	Alp Her	Star	17h 14m 39s	14° 23' 24"	3.5	22h ...	53.7°	149.2°	99.2°	
Rasalhague	Alp Oph	Star	17h 34m 56s	12° 33' 36"	2.1	22h ...	50.1°	143.0°	103.8°	
Thuban	Alp Dra	Star	14h 04m 24s	64° 22' 33"	3.6	01h ...	66.8°	327.3°	83.2°	
Unukalhai	Alp Ser	Star	15h 44m 16s	06° 25' 31"	2.7	00h ...	48.8°	186.6°	75.9°	
Vega	Alp Lyr	Star	18h 36m 56s	38° 46' 55"	0.1	21h ...	60.8°	92.7°	117.7°	
Vindemiatrix	Eps Vir	Star	13h 02m 11s	10° 57' 33"	2.9	02h ...	37.6°	241.0°	39.7°	
Yed Prior	Del Oph	Star	16h 14m 21s	-03° 41' 40"	2.7	23h ...	38.8°	175.9°	81.4°	
HIP73		Star	00h 00m 53s	66° 50' 53"	6.9	16h ...	33.1°	23.9°	121.9°	
HIP128		Star	00h 01m 39s	73° 36' 43"	6.5	15h ...	37.7°	17.9°	115.5°	
HIP207		Star	00h 02m 36s	66° 05' 56"	5.9	15h ...	32.4°	24.4°	122.5°	
HIP274		Star	00h 03m 26s	63° 38' 26"	6.2	15h ...	30.6°	26.3°	124.7°	

Figure 23: Observatory: Catalog tab.




Figure 24: Observatory: All Sky tab.




Figure 25: Observatory: Zoom tab.




Figure 26: **Observatory:** Focus tab with sub-panel *Autofocus Settings* (accessed by the *Options* button).




Figure 27: **Observatory:** Focus tab with sub-panel *Refinement Exposure* (accessed by the *Exposure* button).




Figure 28: **Camera control:** Setup tab with *Setup Camera* sup-panel.




Figure 29: **Camera control:** Setup tab with *Cooler* sup-panel.




Figure 30: **Camera control:** Setup tab with *Setup Filter* sup-panel.




Figure 31: **Camera control:** Expose tab.




Figure 32: **Camera control:** Sub-panel *Autosave Setup* accessible by the Expose tab (see Fig. 31).

6 Spectroscopy with the DADOS spectrograph

6.1 General comments

Operation of the DADOS spectrograph has been tested successfully (see Appendix. B) but the observing procedure and the calibration strategy still need to be improved and established respectively.

Next steps:

- determine approximate correlation between the μm screw positions and observed spectral region
- determine best focus positions for the Slitviewer-Webcam and the spectrum on the ST-1603 camera.
Note: use the calibration lamp and solar spectrum
- Improve and mark alignment of the Webcam with respect to the slits. In principle all slits fit onto the field-of-view of the webcam (see Fig. 39).
Note: align the slits horizontally as in Fig. 39 for N-S and E-W movement of the target.
- align the spectrograph assembly and the ST-1603 camera to the telescope N-S / E-W movement of the target.
- observe different spectral standard stars at different pointing directions and test tracking quality for target on the slit.
- investigate manual tracking options.
Note: for very small pointing corrections of a few arcsec the mount first moves away a greater distance and then slowly comes back to the corrected position (compensation of gear-wheel errors). This is not practical for manual tracking corrections but it might be possible to switch-off or relax this setup for the mount.
- check whether it is more practical to use an external Webcam software (instead of MaxImDL) to simultaneously have the webcam and the telescope controls open
- establish full calibration procedure
Note: see DADOS tutorial by Bernd Koch from Baader planetarium (<http://www.baader-planetarium.de/dados/download/tutorial-dados-d.pdf>)
- check impact of Barlow-Lens, i.e. check vignetting with/without lens.
Note: the spectrograph is optimized for F/10 but the telescope is F/6.8. The barlow lens should correct for that.

6.2 Spectrograph assembly and setup

See Figs. 33-38 for an overview of the packing of all the components and their assembly for observation. See Fig. 41 for the calibration assembly using the neon calibration lamp.

The spectrograph consists of the following sub-parts:

- **DADOS spektrograph:**
 - 3 reflection gratings (200 L/mm, 900 L/mm, 1200 L/mm)
 - 2 eyepieces + focus lens for using eyepieces as slit viewer - different adapter pieces for telescope and camera
- **Webcam** slit viewer + USB Extender-cable
- **ST-1603 CCD camera**
- **1.5x Barlow-lens** to extend telescope focal ratio of F/6.8 to \sim F/10
- **Neon calibration lamp**
- telescope adapter piece

For attachment to the telescope use the adapter piece as shown by Fig. 37. For the **Webcam** use the USB-Extender cable and connect it to the remote telescope computer. For the **ST-1603** use the same USB-cable as for the *STX-16803* which is already laid through the GM4000 mount and connected to the remote computer. Connect the **ST-1603** power adapter to the power outlet box belonging to the main power switch (see Fig. 5).

IMPORTANT: Always make sure that the webcam USB cable and the ST-1603 power cable cannot get stuck when moving the telescope!

IMPORTANT: Make sure that the Barlow-lens is tightly screwed into the DADOS spectrograph and that the complete spectrograph assembly is securely attached to the telescope and all the attachment screws are tight!

IMPORTANT: see DADOS manual for exchanging the gratings. Only exchange the gratings in a clean environment with electrostatic shielding! Use the optical gloves and never touch the gratings by hand!

6.3 Alignment

6.3.1 Spatial orientation

Attach the **ST-1603** and the **Webcam** with orientations as indicated by Fig. 37. This makes sure that all three slits are visible and horizontally oriented in the webcam (see Fig. 39). When attaching the **complete spectrograph assembly** to the telescope make sure that it is oriented horizontal or vertical with respect to N-S / E-W. This makes it easier to move the target onto the slit using the telescope N-S/E-W keys.

Important: Before attaching the spectrograph assembly to the telescope make sure that the focal alignment is done (using the calibration lamp, see next Section)!

6.3.2 Focal alignment

The focal alignment needs to be done off the telescope using the **slit illumination** and the **Neon calibration lamp**:

A) slits on Webcam:

- turn on the **slit illumination** (see Fig. 34).
- **focus the slits** onto the webcam by moving the webcam in/out.
i.e. the slit-image should look similar to Fig. 39 or better
- tightly fix the webcam position
- turn off the slit illumination

B) spectrum on ST-1603:

- Attach the **calibration lamp** to the spectrograph assembly according to Fig. 41.
i.e. replace the Barlow-lens by the 2" receptacle to attach the lamp
- loosen the **Focuser fixation screw** (see Fig. 34)
- start continuous imaging by the ST-1603 camera
- **focus the Neon spectrum** onto the ST-1603 camera by using the camera Focuser of the DADOS spectrograph. (see Figs. 34, 43)
- tighten the Focuser fixation screw
- detach the calibration lamp and re-attach the Barlow-lens

Important: Also make sure that the spatial orientation of the webcam, the camera, and the complete spectrograph assembly is correct (see previous Section)!

6.4 Operation

All sub-components required for the operation of the Spectrograph (i.e. the telescope, the ST-1603 CCD camera, and the Webcam) are controlled by MaxImDL. See Section 5 for operational instructions for the telescope and the CCD camera. In particular note how to setup and connect the ST-1603:

Important: When using the ST-1603 camera set *Filter or Controlling Camera Model* to *No Filter*!

Important: for the ST-1603 *Subframe* does not work and needs to be de-activated!

The webcam is also connected by the **Observatory** window (e.g. see Fig. 20) and operated by the **Webcam** tab (see Fig. 42).

General operation procedure:

- slew telescope to target
- if required: **focus target** onto slit focal plane imaged by ST-1603 by using the **MaxImDL Focus** tab.
 - when using telescope adapter piece as in Fig. 37 focus position $\sim 24913 \mu\text{m}$ (**TBC**)
- **center target onto desired slit** by using the N-S/E-W buttons of the telescope control (see Figs. 21, 42:
 - select step size 10 Sec. or less
 - if required turn On/Off slit illumination
- turn Off slit illumination
- **acquire spectrum** with ST-1603 (see Fig. 44)
- change spectral range by using the grating μm screw (see Fig. 40 and comments therein)
- check target-slit alignment
- **acquire spectrum** with ST-1603
- change spectral range by using the grating μm screw
- check target-slit alignment
- **acquire spectrum** with ST-1603
- ...

6.5 Calibration (TBD)

TBD

see also <http://www.baader-planetarium.de/dados/download/tutorial-dados-d.pdf>

6.6 Figures



Figure 33: Spectrograph assembly: packing.

Top:

(A) DADOS case, (B) Barlow-lens, (C) Webcam, (D) USB Extender-cable, (E) ST-1603 case

Bottom:

(A) Mounting accessories, (B) Eyepieces, (C) Calibration lamp Neon, (D) Barlow-lens, (E) DADOS manual, (F) additional Gratings, (G) DADOS tools, (H) DADOS spectrograph, (J) Webcam, (K) ST-1603 CCD, (L) ST-1603 power cable, (M) ST-1603 accessories




Figure 34: Spectrograph assembly: DADOS spectrograph.
 (A) Slit illumination, (B) towards ST-1603, (C) Focuser for ST-1603 with fixation screw (circle),
 (D) towards telescope, (E) Slit viewer, (F) μm screw for grating position



Figure 35: Spectrograph assembly:
 Connect ST-1603 CCD camera with the DADOS spectrograph.
 (A) ST-1603, (B) adapter piece, (C) DADOS




Figure 36: Spectrograph assembly: Webcam

The webcam has been modified by an adapter piece to attach the slit-viewer focusing lens.

(A) Webcam, (B) adapter piece, (C) Focusing lens




Figure 37: Spectrograph assembly:

Complete assembly with attached Barlow-lens and Webcam ready to be attached to the telescope. Note that this webcam orientation will produce the slit-image as in Fig. 39. The telescope adapter piece as shown is optimized for focusing the target onto the slit.

(A) Barlow-lens, (B) Webcam, (C) Adapter piece for telescope



Figure 38: Spectrograph assembly attached to the telescope.
 Note: in this image the webcam is mounted 180° opposite to the recommended position (see Fig. 37).




Figure 39: Optimized slit orientation on the webcam.




Figure 40: Changing the grating position to change the spectral range:
Loosen the fixation screw (but do not completely remove it!) and use the μm screw to change the illumination angle of the grating. Note that pushing the grating (μm screw goes in) works better than pulling (push fixation screw sideways to support pulling the grating with μm screw):
(A) μm screw, (B) fixation screw




Figure 41: Spectrograph calibration assembly:
Replace the Barlow-lense by the 2" receptacle and attach the Neon calibration lamp. Note: for storage of the Barlow-lens attach it to the barlow-lens adapter piece.
(A) 2" receptacle, (B) Barlow-lens on adapter piece for storage, (C) Neon calibration lamp, (D) telescope adapter piece (not used for calibration)




Figure 42: Webcam operation with *slit illumination* **On** (top) and **Off** (bottom). Settings for the webcam are accessible by the *Options* button.
 Note that in these images the webcam orientation is not in the recommended position (see Figs. 37, 39)




Figure 43: Neon calibration spectrum.




Figure 44: Commissioning image of Vega spectrum.

Note that in this image the webcam and spectrograph assembly orientation are bad. For the recommended orientations see Sect. 6.3 and Figs. 37, 39.

7 Data reduction using MaxImDL (TBD)

TBD

see also video tutorial at http://www.cyanogen.com/maxim_tut.php

A Mounting the SBIG camera

IMPORTANT

Mounting or dismounting the camera requires two people. Do not try it on your own!

DISMOUNTING

Make sure the SBIG power cable and the USB cable are unplugged (see Fig. 48)!
Do not unplug the filter wheel cable and the autoguider cable.

DISMOUNTING

Never dismount the filter wheel which is directly attached to the camera (see Fig. 46)!
For storage attach the camera cover as shown in Fig. 46.

Mounting or dismounting the camera requires two people. One person holds the camera and the other person tightens/loosens the screws. Note that you need to use the US Allen key set (red set, blue set is metric). For using the camera mount the adapter plate that has an edge whereas for visual observations with eye-pieces use the round adapter plate (Allen screw 5/32", see Fig. 45).

- Attach the camera adapter plate to the CDK according to Fig. 45 (with the wedge side facing upwards, use the Allen 1/32" key).
- Loosen all four 5/64" Allen screws at the camera flange (Fig. 46).
- Mount the camera according to Fig. 47 with the autoguider head on the same side as the wedge of the adapter plate. Align the camera parallel to the wedge. Make sure the camera flange tightly fits to the camera adapter plate and then tighten all four 5/64" Allen screws.
- Plugin the SBIG power cable and the USB cable (Fig. 48)




Figure 45: Adapter plates for the camera (*top*) and for the eye-pieces (*bottom*). Note the wedge side facing upwards (away from the CDK control panel). Use the 5/32" Allen key to mount the plates.



Figure 46: SBIG camera with attached filter wheel. For storage attach the cover plate (*top*). For mounting the camera loosen all four 5/64" Allen screws at the camera flange.




Figure 47: Mounting the camera. Note the autoguider head (A) is at the same side as the wedge of the adapter plate. Make sure the camera flange tightly fits to the camera adapter plate and then tighten all four 5/64" Allen screws.




Figure 48: Attach the camera power cable and the USB cable.

Note: in this image instead of the USB cable the Ethernet cable is used. However, only use the USB-cable!

B First light of DADOS spectrograph





Figure 49: Uncalibrated Vega spectrum using the 900 L/mm grating.

The different spectral regions are not calibrated but simply scaled to match each other. The plate scale is ~ 0.096 nm/pixel which was determined using the solar spectrum and the Neon calibration spectrum.



Calibration by Bernd Koch

Figure 50: As comparison to our uncalibrated first light spectrum (Fig. 49):

Flux calibrated Vega spectrum by **Bernd Koch: DADOS Tutorial Baader Planetarium** (<http://www.baader-planetarium.de/dados/download/tutorial-dados-d.pdf>).